



### **Pelican 1660 Carrying Case Specifications**

Logos Imaging provides a custom Pelican 1660 case with an integrated power box and custom foam to hold the scanner, computer, image plates, carousel, and erasing light.

Exterior Case Dimensions: 31.5" L x 22.9" W x 19.0" H  
Case Weight: 45 lbs

Pelican Shipping Cases are Airtight, Watertight up to 30 feet, and they have a Lifetime Warranty. An incredibly lightweight space-age resin, ultra high impact structural copolymer polypropylene, provides four times the strength to weight ratio of conventional injection molding. The cases are unaffected by dents, scratches, or corrosion. They also meet MIL-STD 4150-H stacking (400 lbs.) Drop Test (48"), Immersion (2" at 160°F). Pelican Cases meet and far exceed the highest standards of industrial, airline, military and commercial applications.



### **Custom Soft-Sided Carrying Case**

Logos Imaging also offers a custom soft-sided carrying case option for users that require a lighter weight carrying case with a smaller footprint. The soft-sided case is ideal for use with the T110 scanner with its integrated erasing light and has a custom foam insert to hold the scanner, computer, image plates, and carousel.

Exterior Case Dimensions: 23.5" L x 18.0" W x 18.5" H  
Case Weight: 10.5 lbs